

- ○ QC Wireless SDK Instructions for Use
 - ■ ■ update note:
 - 1.introduce
 - ■ 1.1The role of the SDK
 - 2. API description
 - 2.0 Access Conditions
 - 2.1 Permissions required by the SDK
 - 2.2 Access conditions
 - 2.3.1API
 - Scanning Devices
 - device connection:BleOperateManager.getInstance()
 - 2.3.2 Feature list:
 - Synchronize time, get the list of functions supported by the device
 - bracelet battery
 - Continuous heart rate, blood oxygen, blood pressure switch
 - Set watch sports goals
 - Find equipment
 - Ring photo control
 - Set the Ring to factory reset
 - Wearing Calibration
 - 2.3.3 Data synchronization:
 - Synchronize steps, distance, kcal for the day
 - Synchronized step data details
 - Sync Sleep Data Details
 - Synchronize the details of new sleep data and return according to SetTimeRsp
 - Sedentary data synchronization
 - Sync heart rate data
 - Synchronized blood pressure function
 - Synchronized blood oxygen function
 - Synchronized pressure function
 - Synchronized hrv function

- Synchronized skin temperature function
- Synchronous training records
- 2.3.7 Manual measurement
- 2.3.8 Touch and gestures
- 2.3.9 Changes in Ring measurement data are proactively reported to the APP
- 2.3.10 APP opens exercise type
- Muslim Data Synchronization
-

QC Wireless SDK Instructions for Use

1. Author: James
2. Shenzhen QC.wireless Technology Co., Ltd.
3. Version: 1.0.0

update note:

1. (2021/07/06) scan, connect, measure commands
2. (2021/07/20) Add setting command
3. (2021/07/21) Increase step count, heart rate, sleep data sync
4. (2022/02/28) Add new sleep algorithm
5. (2023/03/10) Add message switch synchronization, body temperature, user information, watch manual measurement of heart rate, blood pressure results
6. (2023/03/16) Add bt connection and contact person
7. (2024/02/23) Add manual pressure, APP movement, pressure synchronization, pressure setting

1.introduce

1.1The role of the SDK

Provide partner companies with the Android Bluetooth SDK for use with Green Orange wireless devices that provide basic and advanced functionality for a major watch or other device. This document is intended to explain the usage context, functionality, etc. of the API. Intended Audience and Reading Recommendations The intended audience and reader recommendations in this article are shown in Annex 1.

Reader	Role
Software Architecture Engineer	Architecture Analysis and Technical Guidance
Android development engineer	Have a certain android development ability, understand Ble related development technology

2. API description

2.0 Access Conditions

Android 5.0 or above, Bluetooth 4.0 or above.

2.1 Permissions required by the SDK

```

//network permissions
<uses-permission android:name="android.permission.INTERNET" />
//Bluetooth related permissions
<uses-permission android:name="android.permission.BLUETOOTH" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADMIN" />
//Storage related permissions
<uses-permission
android:name="android.permission.WRITE_EXTERNAL_STORAGE" />
<uses-permission
android:name="android.permission.READ_EXTERNAL_STORAGE" />
//Location permission
<uses-permission
android:name="android.permission.ACCESS_COARSE_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_BACKGROUND_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_FINE_LOCATION" />

If it is android 12 or higher system
<uses-permission
android:name="android.permission.BLUETOOTH_CONNECT" />
<uses-permission
android:name="android.permission.BLUETOOTH_SCAN" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADVERTISE" />

```

2.2 Access conditions

- Green Orange Wireless Wearables
- Green Orange Wireless SDK and documentation

2.3.1API

Scanning Devices

```

//start scan
BleScannerHelper.getInstance().scanDevice(final Context
context, UUID mUuid, final ScanWrapperCallback scanCallBack);
//stop scan
BleScannerHelper.getInstance().stopScan(Context context)
//Specified device scan
BleScannerHelper.getInstance().scanTheDevice(final Context
context, final String macAddress, final OnTheScanResult scanResult)

```

device connection:BleOperateManager.getInstance()

```

//direct connection

BleOperateManager.getInstance().connectDirectly(smartWatch.deviceAd
dress)
//scan connections

BleOperateManager.getInstance().connectWithScan(smartWatch.deviceAd
dress)
//disconnect
BleOperateManager.getInstance().unBindDevice()
//reconnect
BleOperateManager.getInstance().setNeedConnect(boolean
needConnect)
//Called when bluetooth is turned off
BleOperateManager.getInstance().setBluetoothTurnOff(false)
BleOperateManager.getInstance().disconnect()
//Turn on the system bluetooth monitor
BleOperateManager.getInstance().setBluetoothTurnOff(true)

```

2.3.2 Feature list:

Synchronize time, get the list of functions supported by the device

```

//set time
CommandHandle.getInstance().executeReqCmd(SetTimeReq(0),
ICommandResponse<SetTimeRsp>() {})
//Callback Description
public class SetTimeRsp extends BaseRspCmd {

```

```

//support body temperature
public boolean mSupportTemperature;
//watch face
public boolean mSupportPlate;
//Support the menstrual cycle
public boolean mSupportMenstruation;
//Support custom watch faces
public boolean mSupportCustomWallpaper;
//Support blood oxygen
public boolean mSupportBloodOxygen;
//blood pressure support
public boolean mSupportBloodPressure;
//Support fatigue
public boolean mSupportFeature;
//Support one-key detection
public boolean mSupportOneKeyCheck;
//weather support
public boolean mSupportWeather;
//Support for new sleep protocol
public boolean mNewSleepProtocol;
//Supports up to 6 or 3 dials
public int mMaxWatchFace;
// hrv support
public boolean mSupportHrv;
}

```

Functional support

```

CommandHandle.getInstance()
    .executeReqCmd(
        DeviceSupportReq.getReadInstance(),
        ICommandResponse<DeviceSupportFunctionRsp> {

```

```

//support touch
public boolean supportTouch;
//support muslim
public boolean supportMoslin;
//support ble pair
public boolean supportAPPRvision;
//support Heart rate calibration
public boolean supportBlePair;
//support gesture
public boolean supportGesture;
//support music
public boolean supportRingMusic;
//support video
public boolean supportRingVideo;
//support ebook
public boolean supportRingEbook;

```

```

//support camera
public boolean supportRingCamera;
//support phone call
public boolean supportRingPhoneCall;
//support game
public boolean supportRingGame;
}

```

bracelet battery

```

        CommandHandle.getInstance().executeReqCmd(new
SimpleKeyReq(Constants.CMD_GET_DEVICE_ELECTRICITY_VALUE), new
ICommandResponse<BatteryRsp>() {

            @Override
            public void onDataResponse(BatteryRsp resultEntity) {
                if (resultEntity.getStatus() ==
BaseRspCmd.RESULT_OK) {
                    //Get battery successfully
                }
            }

        });

//Callback Description
public class BatteryRsp extends BaseRspCmd {
//battery [0-100]
private int batteryValue;
}

```

Continuous heart rate, blood oxygen, blood pressure switch

```

//Read continuous heart rate settings
CommandHandle.getInstance()
    .executeReqCmd(HeartRateSettingReq.getReadInstance(),
        ICommandResponse<HeartRateSettingRsp> {
            //it.isEnabled        switch
            //it.heartInterval    heart rate measurement interval

        })

//Read Continuous SpO2 settings
CommandHandle.getInstance()

```

```

        .executeReqCmd(BloodOxygenSettingReq.getReadInstance(),
            ICommandResponse<BloodOxygenSettingRsp> {
                })
        //Read continuous blood pressure settings
        CommandHandle.getInstance()
            .executeReqCmd(BpSettingReq.getReadInstance(),
                ICommandResponse<BpSettingRsp> {
                    })
        //Read pressure setting
        CommandHandle.getInstance()
            .executeReqCmd(PressureSettingReq.getReadInstance(),
                ICommandResponse<PressureSettingRsp>() {
                    switch: it.isEnabled
                })

        //Write continuous heart rate switch isEnabled: true on,
        false: off hrInterval: 10, 15, 20, 30, 60

        CommandHandle.getInstance().executeReqCmd(
            HeartRateSettingReq.getWriteInstance(true,hrInterval),
            ICommandResponse<HeartRateSettingRsp> {
                })
        //Write continuous blood oxygen switch isEnabled: true on,
        false: off
        CommandHandle.getInstance().executeReqCmd(
            BloodOxygenSettingReq.getWriteInstance(boolean isEnabled),
            ICommandResponse<BloodOxygenSettingRsp> {
                })
        //write blood pressure switch

        CommandHandle.getInstance().executeReqCmd(BpSettingReq.getWriteInst
            ance(boolean isEnabled, StartEndTimeEntity startEndTimeEntity, int
            multiple), ICommandResponse<BpSettingRsp> {
                })
        BpSettingRsp,Parameter Description
        isEnabled: true on false off
        StartEndTimeEntity The parameter description is the same as
        above
        multiple default 60

        //Write pressure setting switch,switch isEnabled: true on,
        false: off
        CommandHandle.getInstance().executeReqCmd(

```



```

CommandHandle.getInstance().executeReqCmd(
    PressureSettingReq.getWriteInstance(isEnable),
    ICommandResponse<PressureSettingRsp> {

    })

    //hrv readme
    CommandHandle.getInstance()
        .executeReqCmd(HrvSettingReq.getReadInstance(),
            ICommandResponse<HRVSettingRsp>() {

            })

    //hrv write
    CommandHandle.getInstance().executeReqCmd(
        HrvSettingReq(true),
        ICommandResponse<HRVSettingRsp> {

        })

```

Set watch sports goals

```

CommandHandle.getInstance().executeReqCmd(
    TargetSettingReq.getWriteInstance(
        final int step, final int calorie, final int
distance, final int sportMinute, final int sleepMinute
    ), null
)

```

Parameter Description

step: step target

calorie: Calorie target target, unit card, write kcal to a*1000

distance: Distance to target, in meters

sportMinute: Exercise duration target, in minutes

sleepMinute: Sleep duration target, in minutes

Find equipment

```

CommandHandle.getInstance().executeReqCmd(FindDeviceReq(), null)

```

Ring photo control

```
//The bracelet enters the camera interface

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_INTRO_CAMARA_UI), null)
//The wristband is controlled by the bright screen on the camera interface. The APP will send the bright screen command to keep the watch bright. It is recommended to send it every 2 seconds.
CommandHandle.getInstance().executeReqCmd(

CameraReq(CameraReq.ACTION_KEEP_SCREEN_ON),
            null
        )
//Bracelet click to take a photo event monitoring

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_TAKING_PICTURE,new ICommandResponse<CameraNotifyRsp>(){

    @Override
    public void onDataResponse(CameraNotifyRsp
resultEntity) {

        }
    });
resultEntity.getAction()
Parameter Description
//The watch exits the camera interface
CameraNotifyRsp.ACTION_FINISH
//The watch clicked on the photo event
CameraNotifyRsp.ACTION_TAKE_PHOTO

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_FINISH), null)
```

Set the Ring to factory reset

```
CommandHandle.getInstance().executeReqCmd(RestoreKeyReq(Constants.CMD_RE_STORE),null)
```

- message push

```
//The watch message push switch is fully turned on, and the APP  
should be actively opened  
    CommandHandle.getInstance().executeReqCmd(  
        SetANCSReq(), null  
    )  
//Send message push to watch  
    MessPushUtil.pushMsg(type,message:String)  
  
PushMsgUintReq parameter description  
    type:  
        0x00: Call reminder 0x01: SMS reminder 0x02: QQ reminder  
0x03: WeChat reminder,  
        0x04: incoming call to answer or hang up 0x05: Facebook  
message reminder 0x06: WhatsApp message reminder  
        0x07: Twitter message reminder 0x08: Skype message reminder  
0x09: Line message reminder 0x0a: LinkedIn  
        0x0b: Instagram 0x0c: TIM message 0x0d: Snapchat  
        0x0e: others other types of notifications
```

Wearing Calibration

```

start :enable:true
end :enable:false
    BleOperateManager.getInstance().ringCalibration(false
        ) {
            //2:Measuring 1:success 3:fail
            it.success
        }

    ...

#### firmware version number, hardware version number
```java
 //hardware information

 CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
 .getReadHwRequest());
 //firmware information

 CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
 .getReadFmRequest());

 Receiving implements this QCBluetoothCallbackCloneReceiver
 refer to demo MyBluetoothReceiver
 Judging UUID in the callback onCharacteristicRead

 override fun onCharacteristicRead(uuid: String?, data:
 ByteArray?) {
 if (uuid != null && data != null) {
 val version = String(data, Charsets.UTF_8)
 when(uuid){
 Constants.CHAR_FIRMWARE_REVISION.toString() -> {
 //Firmware version number version
 }
 Constants.CHAR_HW_REVISION.toString() -> {
 //hardware version number version number version
 }
 }
 }
 }
}

```

### 2.3.3 Data synchronization:

Synchronize steps, distance, kcal for the day

```
CommandHandle.getInstance().executeReqCmd(
 SimpleKeyReq(Constants.CMD_GET_STEP_TODAY),
 ICommandResponse<TodaySportDataRsp> {})
```

TodaySportDataRsp parameter description

```
// days ago
private int daysAgo;
// date: year
private int year;
// date: month
private int month;
// date: day
private int day;
// total steps
private int totalSteps;
// running steps/aerobic steps
private int runningSteps;
// calorie value
private int calorie;
// walking distance
private int walkDistance;
// Movement time, in seconds
private int sportDuration;
// sleep time in seconds
private int sleepDuration;
```

## Synchronized step data details

```

 dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
supports synchronization for up to 7 days
 CommandHandle.getInstance().executeReqCmd(
 ReadDetailSportDataReq(dayOffset, 0, 95),
 ICommandResponse<ReadDetailSportDataRsp> {

 })
BleStepDetails parameter description
 //year
 private int year;
 //moon
 private int month;
 //day
 private int day;
 //15 minutes a point, the total number of points in a day is
 96 points, [0, 95], used to calculate the details of the number of
 steps per hour
 private int timeIndex=0;
 // calorie unit card
 private int calorie=0;
 //Step count
 private int walkSteps=0;
 //distance in meters
 private int distance=0;
 // keep for now
 private int runSteps=0;

```

## Sync Sleep Data Details

```

 //deviceAddress Device mac address
 //dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
 supports synchronization for up to 7 days
 //ISleepCallback callback

SleepAnalyzerUtils.getInstance().syncSleepReturnSleepDisplay(device
Address,dayOffset, ISleepCallback {
 //callback: SleepDisplay
})

SleepDisplay parameter description
 // total sleep time
 private int totalSleepDuration;
 // total time of deep sleep
 private int deepDuration;
 // total time of light sleep
 private int shallowDuration;
 // Go to sleep timestamp in seconds
 private int sleepTime;
 // wake up timestamp in seconds
 private int wakeTime;
 // A set of sleep data
 private List<SleepDataBean> list;
 private String address;

SleepDataBean parameter description
 sleepStart start timestamp of a sleep in seconds
 sleepEnd The end timestamp of a sleep in seconds
 type sleep type 2: light sleep 1: deep sleep 3: awake

```

**Synchronize the details of new sleep data and return according to SetTimeRsp**

```

 //offset 0 today 1 yesterday
 public void syncSleepList(int offset, final
ILargeDataSleepResponse response)
 SleepNewProtoResp parameter description
 //sleep start time
 private int st;
 //sleep end time
 private int et;
 //sleep collection
 private List<DetailBean> list;
 DetailBean parameter description
 // duration of a sleep type
 private int d;
 //type of sleep 2: light sleep 3: deep sleep 5: awake
 private int t;

 //support lunch
 public void syncSleepList(int offset, final
ILargeDataSleepResponse response,final
ILargeDataLaunchSleepResponse lunchSleepResponse){

 }
 SleepNewProtoResp
 lunch start time
 private int lunchSt;
 lunch end time
 private int lunchEt;
 private boolean lunchBreak;
 private List<DetailBean> list;

```

## Sedentary data synchronization



```

 //offset 0 today 1 yesterday
 public void syncLongSitList(int offset, final
ILargeSettingForLongDataResponse response)
 LongSitResp parameter description
 // offset 0 today 1 yesterday
 private int index;
 //Sedentary collection
 private List<DetailBean> list;
 DetailBean parameter description
 // duration of a sedentary type
 private int d;
 //0 static (less than 30 steps in 1 minute), 1 trigger
sedentary, 2 movement (more than 30 steps)
 private int t;

```

## Sync heart rate data

```

nowTime current time zone * 3600 + unix second value of current
time
 Sync yesterday: nowTime-(24*3600)*1,
 Sync the day before yesterday: nowTime-(24*3600)*2
 Data can be synchronized for up to three days
 val time = (getTimeZone() * 3600).toInt()
 val nowTime = date.unixTimestamp + time

 CommandHandle.getInstance().executeReqCmd(
 ReadHeartRateReq(nowTime),
 ICommandResponse<ReadHeartRateRsp> {

 })

ReadHeartRateRsp parameter description
 //nothing yet
 private int size = 0;
 //nothing yet
 private int index = 0;
 // unix second value of heart rate data
 private int mUtcTime;
 //The heart rate data array is one point every 5 minutes, the
data subscript *5 is equal to the number of minutes of the day
 private byte[] mHeartRateArray;
 private boolean endFlag = false;

```

## Synchronized blood pressure function

```
//Synchronized automatic blood pressure, measured once an hour
CommandHandle.getInstance()

.executeReqCmd(SimpleKeyReq(Constants.CMD_BP_TIMING_MONITOR_DATA),
ICommandResponse<BpDataRsp> {})
 BpDataEntity parameter description
 //year
 private int year;
 //month
 private int month;
 //day
 private int day;
 private int timeDelay;
 private ArrayList<BpValue> bpValues;

 BpValue parameter description
 //The minute of the day, usually the whole hour
 int timeMinute;
 //measured heart rate value
 int value;

 Get the blood pressure value calculated from the measured value
 //The heart rate value returned by the hr callback, age is the
 age of the user
 val sbp= CalcBloodPressureByHeart.cal_sbp(hr, age) (systolic
 blood pressure)
 //sbp heart rate calculated value
 val dbp=CalcBloodPressureByHeart.cal_dbp(sbp) (diastolic
 pressure)

 //Confirm blood pressure synchronization, call after receiving
 the callback, the watch will delete the records that have been
 synchronized

CommandHandle.getInstance().executeReqCmd(BpReadConformReq(true),nu
ll)

 //synchronize manual blood pressure
 CommandHandle.getInstance()
 .executeReqCmd(ReadPressureReq(0),
ICommandResponse<ReadBlePressureRsp> {})

 ReadBlePressureRsp.getValueList() parameter description
 BpPressure parameter description
 //time seconds value
 public long time;
```

```

public long time;
//(Diastolic pressure)

public int dbp;
//(systolic blood pressure)
public int sbp;

```

## Synchronized blood oxygen function

```

 LargeDataHandler.getInstance().syncBloodOxygenWithCallback(new
IBloodOxygenCallback() {
 @Override
 public void readBloodOxygen(List<BloodOxygenEntity>
data) {

 }
 });

 BloodOxygenEntity parameter description
 //data date
 private String dateStr;
 //Data minimum value array, one data per hour, a total of
24
 private List<Integer> minArray;
 //Data maximum value array, one data per hour, a total of
24
 private List<Integer> maxArray;
 //Data value at 0:00 on a certain day
 private long unix_time;

```

## Synchronized pressure function

```

 //offset 0-6
 //Synchronization pressure 7 days 0 only synchronizes today 1
 yesterday 2 synchronizes the day before yesterday supports up
 to 7 days
 CommandHandle.getInstance()
 .executeReqCmd(PressureReq(offset),
 ICommandResponse<PressureRsp> {

 })

 //PressureRsp Parameter Description
 //Pressure data, the return value divided by 10 is the value
 displayed by the APP, and one data is generated every half hour.
 private byte[] pressureArray;
 //Stress test interval
 private int range=30;
 //Pressure data date
 private DateUtil today;

```

## Synchronized hrv function

```

 //offset 0-6
 //Synchronization hrv 7 days 0 only synchronizes today 1
 yesterday 2 synchronizes the day before yesterday supports up
 to 7 days
 CommandHandle.getInstance()
 .executeReqCmd(HRVReq(0),
 ICommandResponse<HRVRsp> {

 })

 //HRVRsp Parameter Description
 //hrv data, the return value divided by 10 is the value
 displayed by the APP, and one data is generated every half hour.
 private byte[] pressureArray;
 //hrv test interval
 private int range=30;
 //hrv data date
 private DateUtil today;

```

## Synchronized skin temperature function

```

fun registerTempCallback(){
 //Register a callback
 FileHandle.getInstance().clearCallback()
 FileHandle.getInstance().registerCallback(Callback())
 FileHandle.getInstance().initRegister()
}
//Synchronized automatic skin temperature
fun syncAutoTemperature(){
 ///Synchronize automatic body temperature for 3 days 0
 Synchronize only today 1 Synchronize today and yesterday 2
 Synchronize today, yesterday and the day before yesterday
 Support up to 7 days
 FileHandle.getInstance().startObtainTemperatureSeries(2)
}
//Synchronous manual skin temperature
fun syncManual(){
 //Synchronize automatic body temperature for 3 days 0
 Synchronize only today 1 Synchronize today and yesterday 2
 Synchronize today, yesterday and the day before yesterday
 Support up to 7 days
 FileHandle.getInstance().startObtainTemperatureOnce(0)
}

inner class Callback : SimpleCallback() {
 //Continuous temperature callback
 override fun onUpdateTemperature(data: TemperatureEntity) {

 }
 //Single temperature callback
 override fun onUpdateTemperatureList(array:
MutableList<TemperatureOnceEntity>) {

 }
}

TemperatureEntity Parameter Description
public int mIndex; represents today 1 represents yesterday ...6
public int mTimeSpan; Test interval
public float[] mValues; Temperature array cursor*mTimeSpan =
measurement time of the day, in minutes

TemperatureOnceEntity Parameter Description
public long mTime; Measurement timestamp, seconds
public float mValue; Measurement value

```

## Synchronous training records

```
val syncSport= SportPlusHandle()
syncSport.timeFormat="yyyy-MM-dd HH:mm"
syncSport.syncSportPlus { errorCode, t ->

}
//0 is passed for the first synchronization, and the time of the
last movement after synchronization is passed later.
syncSport.cmdSummary(0)

/**
 * Movement type index
 */
public int mSportType;
/**
 * Start time seconds
 */
public int mStartTime;
/**
 * Movement duration seconds
 */
public int mDuration;
/**
 * Movement mileage meters
 */
public int mDistance;
/**
 * Calories small calories
 */
public float mCalories;
/**
 * Average speed cm/s
 */
public int mSpeedAvg;
/**
 * Maximum speed cm/s
 */
public int mSpeedMax;
/**
 * Average heart rate beats/minute
 */
public int mRateAvg;
/**
 * Minimum heart rate beats/minute
 */
public int mRateMin;
```

```

/**
 * Maximum heart rate beats/minute
 */
public int mRateMax;
/**
 * Average altitude cm
 */
public int mElevation;
/**
 * Cumulative climbing cm
 */
public int mUphill;
/**
 * Cumulative downhill cm
 */
public int mDownhill;
/**
 * Average cadence steps/minute
 */
public int mStepRate;
/**
 * Number of exercises times
 */
public int mSportCount;

public int steps;

public List<SportLocation> mLocations = new ArrayList<>();

SportLocation {

/**
 * Real-time heart rate beats/minute
 */
public int mRateReal;
}

typedef NS_ENUM(NSInteger, OdmSportPlusExerciseModelType) {
OdmSportPlusExerciseModelTypeGpsRun = 1, //gps running
OdmSportPlusExerciseModelTypeGpsBike = 2, //gps cycling
OdmSportPlusExerciseModelTypeGpsWalk = 3, //gps walking
OdmSportPlusExerciseModelTypeIndoorRun = 4, //bracelet walking
OdmSportPlusExerciseModelTypeRopeSkipping = 5, //bracelet rope
skipping
OdmSportPlusExerciseModelTypeSwimming = 6, //bracelet swimming
OdmSportPlusExerciseModelTypeRun = 7, //bracelet running-outdoor
OdmSportPlusExerciseModelTypeWalk = 8, //Bracelet hiking
OdmSportPlusExerciseModelTypeBike = 9, //Bracelet cycling

```

```
OdmSportPlusExerciseModelTypeExercise = 10, //Others
OdmSportPlusExerciseModelTypeSwing = 11, //Bracelet swinging
OdmSportPlusExerciseModelTypeClimb = 20, //Bracelet climbing
OdmSportPlusExerciseModelTypeBadminton = 21, //Bracelet badminton
OdmSportPlusExerciseModelTypeYoga = 22, //Bracelet yoga
OdmSportPlusExerciseModelTypeAerobics = 23, //Bracelet aerobics
OdmSportPlusExerciseModelTypeSpinningBike = 24, //Bracelet spinning
bike
OdmSportPlusExerciseModelTypeKayaking = 25, // Bracelet kayaking
OdmSportPlusExerciseModelTypeEllipticalMachine = 26, // Bracelet
elliptical machine
OdmSportPlusExerciseModelTypeRowingMachine = 27, // Bracelet rowing
machine
OdmSportPlusExerciseModelTypePingpong = 28, // Bracelet table
tennis
OdmSportPlusExerciseModelTypeTennis = 29, // Bracelet tennis
OdmSportPlusExerciseModelTypeGolf = 30, // Bracelet golf
OdmSportPlusExerciseModelTypeBasketball = 31, // Bracelet
basketball
OdmSportPlusExerciseModelTypeFootball = 32, // Bracelet football
OdmSportPlusExerciseModelTypeVolleyball = 33, // Bracelet
volleyball
OdmSportPlusExerciseModelTypeRockClimbing = 34, // Bracelet rock
climbing
OdmSportPlusExerciseModelTypeDance = 35, // Bracelet dance
OdmSportPlusExerciseModelTypeRollerSkating = 36, // Bracelet roller
skating

OdmSportPlusExerciseModelTypeTreadmill = 40, // Running-treadmill
OdmSportPlusExerciseModelTypeIndoorWalking = 41, // Running-indoor
walking
OdmSportPlusExerciseModelTypeTrailRunning = 42, // Running-cross-
country running
OdmSportPlusExerciseModelTypeRaceWalk = 43, // Running-race walking
OdmSportPlusExerciseModelTypePlaygroundRunning = 44, // Running -
playground running
OdmSportPlusExerciseModelTypeFatLossRunning = 45, // Running - fat
loss running

OdmSportPlusExerciseModelTypeOutdoorCycling = 50, // Cycling -
outdoor cycling
OdmSportPlusExerciseModelTypeIndoorCycling = 51, // Cycling -
indoor cycling
OdmSportPlusExerciseModelTypeMountainBiking = 52, // Cycling -
mountain biking
OdmSportPlusExerciseModelTypeBMX = 53, // Cycling - BMX

OdmSportPlusExerciseModelTypeSwimmingPool = 55, // Swimming -
```



```
swimming pool swimming
OdmSportPlusExerciseModelTypeOutdoorSwimming = 56, // Swimming -
outdoor swimming
OdmSportPlusExerciseModelTypeFinSwimming = 57, // Swimming-Fin
Swimming
OdmSportPlusExerciseModelTypeSynchronizedSwimming = 58, //
Swimming-Synchronized Swimming

OdmSportPlusExerciseModelTypeOutdoorHiking = 60, // Outdoor Sports-
Outdoor Hiking
OdmSportPlusExerciseModelTypeOrienteering = 61, // Outdoor Sports-
Orienteering
OdmSportPlusExerciseModelTypeFishing = 62, // Outdoor Sports-
Fishing
OdmSportPlusExerciseModelTypeHunt= 63, // Outdoor Sports-Hunting
OdmSportPlusExerciseModelTypeSkateboard = 64, // Outdoor Sports-
Skateboarding
OdmSportPlusExerciseModelTypeParkour = 65, // Outdoor Sports-
Parkour
OdmSportPlusExerciseModelTypeATV = 66, // Outdoor sports - ATV
OdmSportPlusExerciseModelTypeMotocross = 67, // Outdoor sports -
Motocross
OdmSportPlusExerciseModelTypeRacing = 68, // Outdoor sports -
Racing
OdmSportPlusExerciseModelTypeHandCrank = 69,
// Outdoor sports - hand bike
OdmSportPlusExerciseModelTypeMarathon = 70, // Outdoor sports -
marathon
OdmSportPlusExerciseModelTypeObstacleCourse = 71, // Outdoor sports
- obstacle course

OdmSportPlusExerciseModelTypeStairClimber = 80, // Indoor sports -
stair climbing machine
OdmSportPlusExerciseModelTypeStairStepper = 81, // Indoor sports -
stepper
OdmSportPlusExerciseModelTypeMixedAerobic = 82, // Indoor sports -
mixed aerobics
OdmSportPlusExerciseModelTypeKickboxing = 83, // Indoor sports -
kickboxing
OdmSportPlusExerciseModelTypeCoreTraining = 84, // Indoor sports -
core training
OdmSportPlusExerciseModelTypeCrossTraining = 85, // Indoor sports -
cross training
OdmSportPlusExerciseModelTypeIndoorFitness = 86, // Indoor sports -
indoor fitness
OdmSportPlusExerciseModelTypeGroupGymnastics = 87, // Indoor sports
- group gymnastics
OdmSportPlusExerciseModelTypeStrengthTraining = 88, // Indoor
```

```
sports – strength training
OdmSportPlusExerciseModelTypeGapTraining = 89, // Indoor sports –
gap training
OdmSportPlusExerciseModelTypeFreeTraining = 90, // Indoor sports –
free training
OdmSportPlusExerciseModelTypeFlexibilityTraining = 91, // Indoor
sports – flexibility training
OdmSportPlusExerciseModelTypeGymnastics = 92, // Indoor sports –
gymnastics
OdmSportPlusExerciseModelTypeStretch = 93, // Indoor exercise –
stretching
OdmSportPlusExerciseModelTypePilates = 94, // Indoor exercise –
Pilates
OdmSportPlusExerciseModelTypeHorizontalBar = 95, // Indoor exercise
– horizontal bar
OdmSportPlusExerciseModelTypeParallelBars = 96, // Indoor exercise
– parallel bars
OdmSportPlusExerciseModelTypeBattleRope = 97, // Indoor exercise –
battle rope
OdmSportPlusExerciseModelTypeFitness = 98, // Indoor exercise –
fitness
OdmSportPlusExerciseModelTypeBalanceTraining = 99, // Indoor
exercise – balance training
OdmSportPlusExerciseModelTypeStepTraining = 100, // Indoor exercise
– step training

OdmSportPlusExerciseModelTypeSquareDance = 110, // Dance sports–
square dance
OdmSportPlusExerciseModelTypeBallroomDancing = 111, // Dance
sports–social dance
OdmSportPlusExerciseModelTypeBellyDance = 112, // Dance sports–
belly dance
OdmSportPlusExerciseModelTypeBallet = 113, // Dance sports–ballet
OdmSportPlusExerciseModelTypeStreetDance = 114, // Dance sports–
street dance
OdmSportPlusExerciseModelTypeZumba = 115, // Dance sports–Zumba
OdmSportPlusExerciseModelTypeLatinDance = 116, // Dance sports–
Latin dance
OdmSportPlusExerciseModelTypeLatinJazz = 117, // Dance sports–Jazz
dance
OdmSportPlusExerciseModelTypeHipHopDance = 118, // Dance sports –
hip–hop dance
OdmSportPlusExerciseModelTypePoleDancing = 119, // Dance sports –
pole dance
OdmSportPlusExerciseModelTypeBreakDance = 120, // Dance sports –
break dance
OdmSportPlusExerciseModelTypeFolkDance = 121, // Dance sports –
folk dance
```

```
OdmSportPlusExerciseModelTypeNewDance = 122, // Dance sports -
music dance
OdmSportPlusExerciseModelTypeModernDance = 123, // Dance sports -
modern dance
OdmSportPlusExerciseModelTypeDisco = 124, // Dance sports - disco
OdmSportPlusExerciseModelTypeTapDance = 125, // Dance sports - tap
dance
OdmSportPlusExerciseModelTypeOtherDance = 126, // Dance sports -
other dances

OdmSportPlusExerciseModelTypeBoxing = 130, // Combat sports -
boxing
OdmSportPlusExerciseModelTypeWrestling = 131, // Combat sports -
wrestling
OdmSportPlusExerciseModelTypeMartialArts= 132, // Combat sports -
martial arts
OdmSportPlusExerciseModelTypeTaiChi= 133, // Combat sports - Tai
Chi
OdmSportPlusExerciseModelTypeMuayThai= 134, // Combat sports - Muay
Thai
OdmSportPlusExerciseModelTypeJudo= 135, // Combat sports - Judo
OdmSportPlusExerciseModelTypeTaekwondo= 136, // Combat Sports -
Taekwondo
OdmSportPlusExerciseModelTypeKarate= 137, // Combat Sports - Karate
OdmSportPlusExerciseModelTypeFreeSparring= 138, // Combat Sports -
Free Fighting
OdmSportPlusExerciseModelTypeSwordsmanship= 139, // Combat Sports -
Swordsmanship
OdmSportPlusExerciseModelTypeJujitsu= 140, // Combat Sports - Jiu-
Jitsu
OdmSportPlusExerciseModelTypeFencing= 141, // Combat Sports -
Fencing
OdmSportPlusExerciseModelTypeKendo= 142, // Combat Sports - Kendo

OdmSportPlusExerciseModelTypeBeachFootball=150, //!<150Ball Sports
- Beach Soccer
OdmSportPlusExerciseModelTypeBeachVolleyball=151, //!<151Ball
Sports-Beach Volleyball
OdmSportPlusExerciseModelTypeBaseball=152, //!<152Ball Sports-
Baseball
OdmSportPlusExerciseModelTypeSoftball=153, //!<153Ball Sports-
Softball
OdmSportPlusExerciseModelTypeNewFootball=154, //!<154Ball Sports-
Rugby
OdmSportPlusExerciseModelTypeHockey=155, //!<155Ball Sports-Hockey
OdmSportPlusExerciseModelTypeSquash=156, //!<156Ball Sports-Squash
OdmSportPlusExerciseModelTypeDoorKick=157, //!<157Ball Sports-
Gateball
```

OdmSportPlusExerciseModelTypeCricket=158, //!<158Ball Sports-Cricket  
OdmSportPlusExerciseModelTypeHandball=159, //!<159Ball Sports-Handball  
OdmSportPlusExerciseModelTypeBowling=160, //!<160Ball Sports-Bowling  
OdmSportPlusExerciseModelTypePolo=161, //!<161Ball Sports-Polo  
OdmSportPlusExerciseModelTypeRacquetball=162, //!<162Ball Sports-Racquetball  
OdmSportPlusExerciseModelTypeBilliards=163, //!<163Ball Sports-Table Tennis  
OdmSportPlusExerciseModelTypeTakraw=164, //!<164Ball Sports-Sepak Rattan Ball  
OdmSportPlusExerciseModelTypeDodgeBall=165, //!<165Ball Sports-Dodgeball  
OdmSportPlusExerciseModelTypeWaterPolo=166, //!<166Ball Sports-Water Polo  
OdmSportPlusExerciseModelTypePuck=167, //!<167Ball Sports-Ice Hockey  
OdmSportPlusExerciseModelTypeShuttlecock=168, //!<168Ball Sports-Shuttlecock  
OdmSportPlusExerciseModelTypeIndoorSoccer=169, //!<169Ball Sports-Indoor Football  
OdmSportPlusExerciseModelTypeSandbag=170, //!<170Ball Sports-Sandbag Ball  
OdmSportPlusExerciseModelTypeBocce=171, //!<171Ball Sports-Floor Bocce  
OdmSportPlusExerciseModelTypeJaiBall=172, //!<172Ball Sports-Paiai  
OdmSportPlusExerciseModelTypeFloorBall=173, //!<173Ball Sports-Floorball  
OdmSportPlusExerciseModelTypeAustralianRulesFootball=174, //!<174Ball Sports-Australian Rules Football  
OdmSportPlusExerciseModelTypePickering=175, //!<175Ball Sports-Pickering  
  
OdmSportPlusExerciseModelTypeOutdoorBoating=180, //!<180Water Sports-Outdoor Rowing  
OdmSportPlusExerciseModelTypeSailing=181, //!<181Water Sports-Sailing  
OdmSportPlusExerciseModelTypeDragonBoat=182, //!<182 Water Sports-Dragon Boat  
OdmSportPlusExerciseModelTypeSurf=183, //!<183 Water Sports-Surfing  
OdmSportPlusExerciseModelTypeKitesurfing=184, //!<184 Water Sports-Kite Surfing  
OdmSportPlusExerciseModelTypePaddling=185, //!<185 Water Sports-Paddling  
OdmSportPlusExerciseModelTypePaddleboard=186, //!<186 Water Sports-Paddleboard Surfing

OdmSportPlusExerciseModelTypeIndoorSurfing=187, //!<187 Water Sports-Indoor Surfing  
OdmSportPlusExerciseModelTypeDrifting=188, //!<188 Water Sports-Rafting  
OdmSportPlusExerciseModelTypeSnorkeling=189, //!<189 Water Sports - Snorkeling

OdmSportPlusExerciseModelTypeSkis=190, //!<190 Snow Sports - Double Skiing  
OdmSportPlusExerciseModelTypeSnowboard=191, //!<191 Snow Sports - Singleboarding  
OdmSportPlusExerciseModelTypeAlpineSkiing=192, //!<192 Snow Sports - Alpine Skiing  
OdmSportPlusExerciseModelTypeCrossCountrySkiing=193, //!<193 Snow Sports - Cross-Country Skiing  
OdmSportPlusExerciseModelTypeSkiOrientserlng=194, //!<194 Snow Sports - Orienteering  
OdmSportPlusExerciseModelTypeBiathlon=195, //!<195 Snow Sports - Biathlon  
OdmSportPlusExerciseModelTypeOutdoorSkating=196, //!<196Snow and ice sports-outdoor skating  
OdmSportPlusExerciseModelTypeIndoorSkating=197, //!<197Snow and ice sports-indoor skating  
OdmSportPlusExerciseModelTypeCurling=198, //!<198Snow and ice sports-curling  
OdmSportPlusExerciseModelTypeBobsleigh=199, //!<199Snow and ice sports-bobsleigh  
OdmSportPlusExerciseModelTypeSled=200, //!<200Snow and ice sports-sled  
OdmSportPlusExerciseModelTypeSnowmobile=201, //!<201Snow and ice sports-snowmobile  
OdmSportPlusExerciseModelTypeSnowshoeing=202, //!<202 Ice and Snow Sports-Snowshoe Hiking

OdmSportPlusExerciseModelTypeHulaHoop=210, //!<210 Leisure Sports-Hula Hoop  
OdmSportPlusExerciseModelTypeFrisbee=211, //!<211 Leisure Sports-Frisbee  
OdmSportPlusExerciseModelTypeDarts=212, //!<212 Leisure Sports-Darts  
OdmSportPlusExerciseModelTypeFlyAKite=213, //!<213 Leisure Sports-Kite Flying  
OdmSportPlusExerciseModelTypeTugOfWar=214, //!<214 Leisure Sports-Tug of War  
OdmSportPlusExerciseModelTypeEsports=215, //!<215 Leisure Sports-Esports  
OdmSportPlusExerciseModelTypeStroller=216, //!<216 Leisure Sports-Walking Machine

```

OdmSportPlusExerciseModelTypeNewSwing=217, //!<217 Leisure Sports-
Swing
OdmSportPlusExerciseModelTypeShuffleboard=218, //!<218 Leisure
Sports-Shuffleboard
OdmSportPlusExerciseModelTypeTableSoccer=219, //!<219 Leisure
Sports-Table Football
OdmSportPlusExerciseModelTypeSomatosensoryGame=220, //!<220 Leisure
Sports-Somatics
OdmSportPlusExerciseModelTypeBungeeJumping=221, //!<221 Leisure
Sports-Bungee Jumping
OdmSportPlusExerciseModelTypeParachute=222, //!<222 Leisure Sports-
Skydiving
OdmSportPlusExerciseModelTypeAnusara=223, //!<223 Leisure Sports-
Anusara
OdmSportPlusExerciseModelTypeYinYoga=224, //!<224 Leisure Sports-
Yin Yoga
OdmSportPlusExerciseModelTypePregnancyY
``

```

### 2.3.6 OTA upgrade function:

```

```java
//dfu upgrade instance
val fuHandle= DfuHandle.getInstance()
//initialize callback
dfuHandle.initCallback()
//DFU file verification, path firmware file path
if (dfuHandle.checkFile(path)) {
    dfuHandle.start(dfuOpResult)
}
//dfuOpResult callback description
private val dfuOpResult: DfuHandle.IOpResult = object :
DfuHandle.IOpResult {
    override fun onActionResult(type: Int, errCode: Int) {
        if (errCode == DfuHandle.RSP_OK) {
            when (type) {
                1 -> dfuHandle.init()
                2 -> dfuHandle.sendPacket()
                3 -> dfuHandle.check()
                4 -> {
                    //The upgrade is successful, wait for the
device to restart

                    dfuHandle.endAndRelease()
                }
            }
        } else {
            //Upgrade exception or failure
        }
    }
}
```

```

```

 }

 override fun onProgress(percent: Int) {
 // file upgrade progress
 }
}

```

## 2.3.7 Manual measurement

```

//StartHeartRateRsp parameter description
 private byte type; type 1: heart rate 2: blood
pressure 3: blood oxygen
 private byte errCode; measurement error code 0:
normal 1: measurement failed 2: measurement failed
 private byte value; measurement value: heart rate
or blood oxygen
 private byte sbp; blood pressure sbp
 private byte dbp; blood pressure dbph

 Boolean stop true stop false start

// manual heart rate
 BleOperateManager.getInstance().manualModeHeart(new
ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp
resultEntity) {

 }
 },Boolean stop);
//manual blood pressure
 BleOperateManager.getInstance().manualModeBP(new
ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp
resultEntity) {

 }
 },Boolean stop);
//manual blood oxygen
 BleOperateManager.getInstance().manualModeSpO2(new
ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp

```

```

resultEntity) {

 }
 },Boolean stop);
 //manual pressure
 BleOperateManager.getInstance().manualModePressure(new
 ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp
resultEntity) {

 }
 },Boolean stop);
 //manual hrv
 BleOperateManager.getInstance().manualModeHrv(new
 ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp
resultEntity) {

 }
 },Boolean stop);
 //manual temperature
 BleOperateManager.getInstance().manualTemperature(new
 ICommandResponse<StartHeartRateRsp>() {
 @Override
 public void onDataResponse(StartHeartRateRsp
resultEntity) {
 //The temperature value obtained should be divided by
 10 to get the normal temperature.
 }
 },Boolean stop);
 }
 }
 }
 }
 }
}

```

## 2.3.8 Touch and gestures



```

 read
 touch: true touch false:gestures

CommandHandle.getInstance().executeReqCmd(TouchControlReq.getReadIn
stance(false),
 ICommandResponse<TouchControlResp> {
 appType: 0:close 1:music 2:video 3:muslim 4:ebook
5:camera 6:phone call 7:game 8:heart
 strength: Dynamics
 }

 write:
 appType: 0:close 1:music 2:video 3:muslim 4:ebook 5:camera
6:phone call 7:game 8:heart
 touch: true touch false:gestures
 Strength:1-10

CommandHandle.getInstance().executeReqCmdNoCallback(TouchControlReq
.getWriteInstance(appType,false,currStrength))

```

## 2.3.9 Changes in Ring measurement data are proactively reported to the APP

```

//Add listener

BleOperateManager.getInstance().addOutDeviceListener(ListenerKey.He
art,myDeviceNotifyListener)

ListenerKey Parameter Description
public class ListenerKey {
 public static int Heart=1; Heart
 public static int BloodPressure=2; Blood Pressure
 public static int BloodOxygen=3; Blood Oxygen
 public static int Temperature=5; Temperature
 public static int SportRecord=7; Sport Record
 public static int All=7; All
}

Monitoring instructions
 inner class MyDeviceNotifyListener : DeviceNotifyListener() {
 override fun onDataResponse(resultEntity: DeviceNotifyResp?)
 {
 if (resultEntity!!.status == BaseRspCmd.RESULT_OK) {

```

```

BleOperateManager.getInstance().removeOthersListener()
 when (resultEntity.dataType) {
 1 -> {
 //Watch heart rate test changes
 }
 2 -> {
 //Watch blood pressure test changes
 }
 3 -> {
 //Watch blood oxygen test changes
 }
 4 -> {
 //Changes in watch step counting details
 }
 5 -> {
 //Watch body temperature changes on the day
 }
 7 -> {
 //Generate new exercise records
 }
 0x0c -> {
 //730c5701000000000000000000000000d7
 val charging =
 BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[2]
)
)
 if (charging > 0) {
 //charging
 } else {
 val battery =
 BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[1]
)
)
 //battery power
 }
 }
 0x12 -> {
 //7312 00005200025100003c0000000066
 AwLog.i(Author.HeZhiYuan,
 ByteUtil.bytesToString(resultEntity.loadData))

 val step = BLEDataFormatUtils.bytes2Int(

```

```

 byteArrayOf(
 resultEntity.loadData[1],
 resultEntity.loadData[2],
 resultEntity.loadData[3]
)
)
)

 val calorie = BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[4],
 resultEntity.loadData[5],
 resultEntity.loadData[6]
)
)
)
 val distance =
BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[7],
 resultEntity.loadData[8],
 resultEntity.loadData[9]
)
)
)
 deviceNotification(step, distance, calorie)
 }
 }
 }sl
 }
 }

 //Remove the listener. It must be removed after registration,
 otherwise multiple callbacks will appear.

 BleOperateManager.getInstance().removeNotifyListener(ListenerKey.He
art)
 //Remove all listeners

 BleOperateManager.getInstance().removeNotifyListener(ListenerKey.Al
l)

```

## 2.3.10 APP opens exercise type

```

 // status 1 Start movement 2 Pause 3 Continue 4 End 6 Movement
start timestamp
 //Sport type 4 Walking 5 Jumping rope 7 Running 8 Hiking 9
Cycling 10 Other sports 20 Hiking 21 Badminton
 22 Yoga 23 Aerobics 24 Spinning 25 Kayaking 26 Elliptical

```

machine 27 Rowing machine 28 Table tennis 29 Tennis

30 Golf 31 Basketball 32 Football 33 Volleyball 34 Rock  
climbing 35 Dance 36 Roller skating 60 Outdoor hiking

```
CommandHandle.getInstance().executeReqCmd(
 PhoneSportReq.getSportStatus(
 1, sportType
), gpsResponse
)
```

```
private var gpsResponse: ICommandResponse<AppSportRsp> =
 ICommandResponse<AppSportRsp> { resultEntity ->
 AwLog.i(Author.HeZhiYuan, resultEntity)
 if (resultEntity != null) {
 when (resultEntity.gpsStatus) {
 6 -> {
 //Exercise start time (Unit second)
 }

 2 -> {
 //Exercise pause
 }

 3 -> {
 // //Exercise continues
 }

 4 -> {
 //Exercise end
 }

 0x25 -> {
 //Muslim ring click real-time data
 //732500000013000700000000000000b2
 val count = BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[1],
 resultEntity.loadData[2],
 resultEntity.loadData[3],
 resultEntity.loadData[4],
)
)
 }
 }
 }
 }
```

//Report data during exercise

//Add motion data reporting and monitoring

BleOperateManager.getInstance().addSportDeviceListener(0x78

```

BleOperateManager.getInstance().addSportDeviceListener(0x78,
myDeviceSportNotifyListener)

//Remove sports data reporting monitoring

BleOperateManager.getInstance().removeSportDeviceListener(0x78)
//Listening to inner classes
inner class MyDeviceNotifyListener :
DeviceSportNotifyListener() {
 override fun onDataResponse(resultEntity: DeviceNotifyRsp?)
{
 super.onDataResponse(resultEntity)

 if (resultEntity!!.status ==
BaseRspCmd.RESULT_OK) {
 //Movement duration, unit seconds
 val sportTime = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[2],
 resultEntity.loadData[3]
)
)
 //Exercise real-time heart rate
 val heart = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[4]
)
)

 //The number of steps generated during exercise
will only have data when the exercise type is 4, 7, or 8, otherwise
it will be 0
 val step = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[5],
 resultEntity.loadData[6],
 resultEntity.loadData[7]
)
)
 //The distance generated during exercise will only
have data when the exercise type is 4, 7, or 8, and the others will
be 0.
 val distance = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[8],
 resultEntity.loadData[9],
 resultEntity.loadData[10]
)
)
 //Calories generated during exercise

```

```

 //calories generated during exercise
 val calorie = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[11],
 resultEntity.loadData[12],
 resultEntity.loadData[13]
)
)

 //Error status during exercise
 val status = bytes2Int(
 byteArrayOf(
 resultEntity.loadData[1]
)
)
 val sportType = BLEDataFormatUtils.bytes2Int(
 byteArrayOf(
 resultEntity.loadData[0]
)
)
 if (status == 0x03) {
 //Not wearing
 }
 }
}

/**
 * Convert the byte array to int type, with the high byte of
the array first
 *
 * @param data
 * @return
 */
public static int bytes2Int(byte[] data) {
 int length = data.length;
 int res = 0;
 for (int i = 0; i < length; i++) {
 res |= (data[i] & 0xFF) << (8 * (length - 1 - i));
 }
 return res;
}

```

## Muslim Data Synchronization

```
//dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
supports synchronization for up to 7 days
CommandHandle.getInstance()
 .executeReqCmd(
 MuslimReq(0),
 ICommandResponse<MuslimRsp> {

 })

//MuslimRsp
private int size = 0;
private int index = 0;
private byte[] pressureArray; //Muslim data one data point per
hour
private boolean endFlag = false;
private int range=60;
private DateUtil today;
private int offset=-1;
```

Real-time data reporting, ring single reporting, please refer to  
2.3.9 Type 0x25